

BC-523A

Enriched Boron-loaded Liquid Scintillator

BC-523A liquid scintillator contains boron enriched in ^{10}B to a minimum of 90%. The scintillator is useful in total absorption neutron spectrometry.

A fast neutron will produce a prompt recoil proton pulse with its initial scatterings in the liquid. Neutrons that are sufficiently thermalized within the scintillator are likely to undergo the ^{10}B (n, α) capture. The capture pulse is in delayed coincidence with the prompt pulse and is used to identify neutron events. For neutron energies below 200 keV, the capture time constant is determined solely by the ^{10}B concentration. The average capture time is about 0.5 microseconds.

The scintillator also exhibits excellent pulse shape discrimination properties for neutron-gamma separation. References dealing with the performance of boron-loaded organic scintillators are listed on the back.

BC-523A is supplied sealed in ready-to use glass or metal cells. The liquid is sensitive to moisture, and you must take care to avoid exposing it to air if you perform your own encapsulation.

General Technical Data –

^{10}B Content	4.41%
Density	0.916 g/cc
Refractive Index	1.415
Flash Point	-8°C

Scintillation Properties

Light Output, % Anthracene	65
Decay time, short component, ns	3.7
Wavelength of maximum emission, nm	425
Bulk light attenuation, meters	>4

Atomic Composition –

No. of H Atoms per cc	4.98×10^{22}
No. of C Atoms per cc	2.86×10^{22}
No. of O Atoms per cc	0.811×10^{22}
No. of ^{10}B Atoms per cc	0.243×10^{22}
No. of ^{11}B Atoms per cc	0.027×10^{22}

BC-523A Natural Boron-loaded Liquid Scintillator

Emission Spectrum –



USA

Saint-Gobain Crystals
17900 Great Lakes Parkway
Hiram, OH 44234
Tel: (440) 834-5600
Fax: (440) 834-7680

Europe

Saint-Gobain Crystals
104 Route de Larchant
BP 521
77794 Nemours Cedex, France
Tel: 33 (1) 64 45 10 10
Fax: 33 (1) 64 45 10 01

P.O. Box 3093
3760 DB Soest
The Netherlands
Tel: 31 35 60 29 700
Fax: 31 35 60 29 214


Japan

Saint-Gobain KK, Crystals Division
3-7, Kojimachi, Chiyoda-ku,
Tokyo 102-0083 Japan
Tel: 81 (0) 3 3263 0559
Fax: 81 (0) 3 5212 2196

China

Saint-Gobain (China) Investment Co., Ltd.
15-01 CITIC Building
19 Jianguomenwai Ave.
Beijing 100004 China
Tel: 86 (0) 10 6513 0311
Fax: 86 (0) 10 6512 9843

www.detectors.saint-gobain.com



References:

1. L.W. Bollinger and G.E. Thomas, Rev. Sci. Instrum., 28, 489-496 (1957)
2. Von L. Sutterlin, Atomkernenergie, Vol. 12, NO. 7/8, 287,288 (1967)
3. L.R. Greenwood and N.R. Chellew, Rev. Sci. Instrum., 50 (4), 466-471 (April, 1979)
4. D.M. Drake, et al, Nucl. Instr. & Methods in Phy. Res., A274 576-582 (1986)