

S350C Thermal Power Head

The S350C is a power meter head, designed for power measurements on industrial Excimer lasers up to 60W with low drift and high resolution. The sensor has a large active area with 40mm diameter and a high damage threshold excimer coating.

To perform accurate measurements the S350C has to be zeroed before starting a measurement. Further it is recommended to prevent the sensor from air flow or other thermal disturbances. The head has a built in thermistor to control a overheating of the sensor.

The S310C is compatible with all available Thorlabs power meter consoles. A non-volatile memory in the sensor connector contains sensor information data and the NIST and PTB traceable calibration data.

Technical Specifications


Sensor Model	S350C
Optical Power Range (Continuous)	10mW - 40W
Wavelength Range	0.19 - 1.1µm & 10.6µm
Application	Pulsed, long pulse (Excimers @ 190nm, 258nm, 351nm), Tea, CO ₂ , YAG 3th & 4th)
Detector Type	Thermal Surface Absorber
Coating	Excimer EXC
Detector Size (active area)	Ø40mm
Distance Frontface to Detector	13mm
Power Resolution ¹⁾	1mW
Calibration Uncertainty	+/- 3% @351nm +/- 5% @ 190nm - 1100nm, 10.6µm
Linearity	+/- 1%
Max Intermittent Power (2min max)	60W
Max. Average Power Density	2kW/cm ²
Max. Pulse Energy Density	0.7J/cm ² (1ns pulse) 10J/cm ² (1ms pulse)
Response time with display (0-90%) ²⁾	1 sec.
Cooling	Convection
Head Temperature Measurement	Thermistor 2kΩ
Sensor Dimensions	100 x 100 x 55 mm ³
Connector Cable Length	1.5m
Connector	Sub-D 9p male
Mounting and accessories	N/A
Post	M6, 75mm post included
Weight	1 kg
Console Compatibility	PM100D, PM100A, PM100USB, PM200, PM320E

¹⁾ with PM100D console, acceleration circuit off

²⁾ Other calibration wavelengths on request

³⁾ Spectral calibration - the response values for wavelength correction outside this calibration range, is interpolated from the general absorption curve of the absorber.

Mechanical Drawing



Available Accessories


The S350C is compatible to the Thorlabs imperial and metric post and post-holder series.

Cleaning and Maintenance


There are no serviceable parts in the S350C thermal head. The housing may be cleaned by wiping with a soft damp cloth. The detectors on the thermal heads cannot be cleaned. Gently blow off any debris using compressed air. If any scratches or other signs of damage remain on the sensor area, contact Thorlabs service department for repair or replacement. If you suspect a problem with your S350C please call Thorlabs and an engineer will be happy to assist you.

As long as the sensor has not been exposed to excessive optical power (please pay attention to the maximum ratings in the technical specifications), the calibration should be very stable over long periods of time (well over a year). To keep the accuracy and performance of the S350C, Thorlabs recommends a yearly recalibration, starting one year after purchase.

General Absorption Curves



Pulse Energy Ratings and Damage Thresholds



WEEE

As required by the WEEE (Waste Electrical and Electronic Equipment Directive) of the European Community and the corresponding national laws, Thorlabs offers all end users in the EC the possibility to return “end of life” units without incurring disposal charges.

This offer is valid for Thorlabs electrical and electronic equipment

- sold after August 13th 2005
- marked correspondingly with the crossed out “wheelie bin” logo (see fig. 1)
- sold to a company or institute within the EC
- currently owned by a company or institute within the EC
- still complete, not disassembled and not contaminated

As the WEEE directive applies to self contained operational electrical and electronic products, this “end of life” take back service does not refer to other Thorlabs products, such as

- pure OEM products, that means assemblies to be built into a unit by the user (e. g. OEM laser driver cards)
- components
- mechanics and optics
- left over parts of units disassembled by the user (PCB's, housings etc.).

If you wish to return a Thorlabs unit for waste recovery, please contact Thorlabs or your nearest dealer for further information.


Waste treatment on your own responsibility

If you do not return an “end of life” unit to Thorlabs, you must hand it to a company specialized in waste recovery. Do not dispose of the unit in a litter bin or at a public waste disposal site.

Ecological background

It is well known that WEEE pollutes the environment by releasing toxic products during decomposition. The aim of the European RoHS directive is to reduce the content of toxic substances in electronic products in the future.

The intent of the WEEE directive is to enforce the recycling of WEEE. A controlled recycling of end of live products will thereby avoid negative impacts on the environment.



Crossed out “wheelie bin” symbol