### 387/480 nm BrightLine® dual-band bandpass filter

## Part Number: FF01-387/480-25





#### Semrock, Inc.

3625 Buffalo Road, Suite 6 Rochester, New York 14624

Main Phone: +1 585.594.7050 (worldwide)
Toll Free Phone: 866.736.7625 (866-SEMROCK)
(within US and Canada)

Your filter spectrum may differ slightly from the typical spectrum above, but is certified to meet the optical specifications noted below.



#### 387/480 nm BrightLine® dual-band bandpass filter

Individual multiband fluorescence bandpass filters that utilize Semrock's patented single-substrate construction. These filters have extremely high transmission, steep and well-defined edges, and outstanding blocking between the passbands. All thin-film, hard-coated construction for unsurpassed performance and reliability.

Part Number	Size	Price1	Stock Status
FF01-387/480-25	25 mm x 5.0 mm	\$395	In Stock
FF01-387/480-21.8-D	21.8 mm x 2.0 mm (unmounted)	\$395	2nd Day Ship

Don't see a size you need? Contact us for custom sizing - available in less than a week (sizing fee applies).

1) US domestic pricing only. If you are ordering from outside the US, please contact your nearest regional distributor for the correct list price.

### **Optical Specifications**

Value
Tavg > 80% 381.5 – 392.5 nm
387 nm
11 nm
15.8 nm
Tavg > 90% 466 - 495 nm
480 nm
29 nm
34.3 nm
ODavg > 6 300 - 374.4 nm
OD > 3.5 403.5 nm
ODavg > 6 414 - 452 nm
ODavg > 3.5 459 - 502.5 nm
ODavg > 6 410 - 550 nm
ODavg > 6 550 - 725 nm
ODavg > 4 725 - 850 nm

## **General Filter Specifications**

Specification	Value	
Angle of Incidence	0 ± 5 degrees	
Cone Half-angle	7 degrees	
Optical Damage Rating	Testing has proven to show no signs of degradation when exposed to at least 6.0 W of power from an unfiltered xenon arc lamp over a 25 mm diameter (corresponding to 1.2 W/cm²) for over 500 hrs.	
Effective Index	1.91	

# Physical Filter Specifications (applies to standard sized parts; contact us regarding other sizes)

Specification	Value
Transverse Dimensions (Diameter)	25 mm
Transverse Tolerance (mounted)	+ 0.0 / – 0.1 mm
Filter Thickness (Mounted)	5.0 mm
Filter Thickness Tolerance (Mounted)	± 0.1 mm
Clear Aperture	≥ 21 mm
Scratch-Dig	60-40
Substrate Thickness (unmounted)	2.0 mm
Substrate Thickness Tolerance (unmounted)	± 0.1 mm
Orientation	Arrow on ring indicates preferred direction of propagation of light